



# FAZZINI

Your Health is Our Aim

Fazzini srl  
S.S Padana Superiore, 317  
20055 - Vimodrone-Milano-ITALY  
Tel +39022651521- Fax +390226515221  
info@fazzini.it - www.fazzini.it

## AUTOMATIC LOADING STRETCHER WITH MATTRES

910A/W



### DIMENSION:

191 x 55 x 83 H cm

Load 159 kg



### GENERAL DESCRIPTION:

Made of anodized aluminum  
Mounted on 4 wheels of 125 mm diameter, 2 fixed and 2 rotating self-centering  
Loading height regulation system allows an easy loading on the ambulance  
Adjustable backrest 0° - 75°  
Foldable side rails.  
Supplied with washable mattress  
Locking system for folded position  
Safety lever situated near the unlocking device



### STANDARD ACCESSORIES:

Washable mattress  
Blocking system  
2 belts  
Nylon bumpers for vertical legs  
125mm castors two rotating  
Side protections completely overturnable

### ACCESSORIES ON DEMAND

I.V. stand  
oxygen tank holder

