

USE

For daily catheterisation, for short-term use for introduction through the urethra for urine collection or for therapeutic purposes.

CHARACTERISTICS

Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining

- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter
 - Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections
 - Water sachet attached, for easy activation of the coating
 - Some dimensions have an additional sleeve allowing aseptic handling



STERILE EO



0123

Catalogue No.	Colour	CH	Length	Outer diameter	Note	Packaging
62294	Blue	8	400 mm	2.6 mm	Male	30 pieces
62295	Black	10	400 mm	3.3 mm	Male	30 pieces
62296	White	12	400 mm	4.0 mm	Male - with sleeve	30 pieces
62297	Green	14	400 mm	4.6 mm	Male - with sleeve	30 pieces
62298	Orange	16	400 mm	5.3 mm	Male	30 pieces
62299	Red	18	400 mm	6.0 mm	Male	30 pieces
62300	Blue	8	180 mm	2.6 mm	Female	30 pieces
62301	Black	10	180 mm	3.3 mm	Female	30 pieces
62302	White	12	180 mm	4.0 mm	Female	30 pieces
62303	Green	14	180 mm	4.6 mm	Female	30 pieces
62304	Orange	16	180 mm	5.3 mm	Female	30 pieces
62305	Red	18	180 mm	6.0 mm	Female	30 pieces
62306	Green	6	280 mm	2.0 mm	Pediatric	30 pieces
62307	Blue	8	280 mm	2.6 mm	Pediatric - with sleeve	30 pieces
62308	Black	10	280 mm	3.3 mm	Pediatric - with sleeve	30 pieces

greencath[®]
PLUS



Coated Nelaton
urethral catheter
with water sachet

