

# greencath® PLUS

## Coated Nelaton urethral catheter with water sachet

- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter
- Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections
- Water sachet attached, for easy activation of the coating, containing specially purified, sterile water
- Some dimensions have an additional sleeve allowing aseptic handling

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
62294	■ blue	8	400	2,6	Male	30	600
62295	■ black	10	400	3,3	Male	30	600
62296*	□ white	12	400	4,0	Male	30	600
62297*	■ green	14	400	4,6	Male	30	600
62298	■ orange	16	400	5,3	Male	30	600
62299	■ red	18	400	6,0	Male	30	600
62300	■ blue	8	180	2,6	Female	30	900
62301	■ black	10	180	3,3	Female	30	900
62302	□ white	12	180	4,0	Female	30	900
62303	■ green	14	180	4,6	Female	30	900
62304	■ orange	16	180	5,3	Female	30	900
62305	■ red	18	180	6,0	Female	30	900
62306	■ green	6	280	2,0	Paediatric	30	600
62307*	■ blue	8	280	2,6	Paediatric	30	600
62308*	■ black	10	280	3,3	Paediatric	30	600

\* - with sleeve



POLISHED SIDE EYES



WATER SACHET



HYDROPHILIC COATING

