



greencath®SOFT

Coated soft tip urethral catheter


- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter
- Polished side eyes and tube flexibility additionally reduce friction between the catheter and urethral mucosal lining
- Soft, flexible tip enables easy catheter introduction in case of enlarged prostate and various urethral masses



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
62549	■ blue	8	400	2,6	Male	30	600
62550	■ black	10	400	3,3	Male	30	600
62551	□ white	12	400	4,0	Male	30	600
62552	■ green	14	400	4,6	Male	30	600
62553	■ orange	16	400	5,3	Male	30	600
62554	■ blue	8	180	2,6	Female	30	900
62555	■ black	10	180	3,3	Female	30	900
62556	□ white	12	180	4,0	Female	30	900
62557	■ green	14	180	4,6	Female	30	900
62558	■ orange	16	180	5,3	Female	30	900
62560	■ blue	8	310	2,6	Paediatric	30	600
62561	■ black	10	310	3,3	Paediatric	30	600



POLISHED
SIDE EYES



SOFT
TIP



HYDROPHILIC
COATING