




# greencath®

## Coated Nelaton urethral catheter



### USE

For daily catheterisation, for short-term use for introduction through the urethra for urine collection or for therapeutic purposes.

### CHARACTERISTICS

Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining

- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter •

Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections

STERILE EO



HYDROPHILIC  
COATING

Catalogue No.	Colour	CH	Length	Outer diameter	Note	Packaging
59140	Blue	8	400 mm	2.6 mm	Male	30 pieces
59141	Black	10	400 mm	3.3 mm	Male	30 pieces
59142	White	12	400 mm	4.0 mm	Male	30 pieces
59143	Green	14	400 mm	4.6 mm	Male	30 pieces
59144	Orange	16	400 mm	5.3 mm	Male	30 pieces
59145	Red	18	400 mm	6.0 mm	Male	30 pieces
59146	Blue	8	180 mm	2.6 mm	Female	30 pieces
59147	Black	10	180 mm	3.3 mm	Female	30 pieces
59148	White	12	180 mm	4.0 mm	Female	30 pieces
59149	Green	14	180 mm	4.6 mm	Female	30 pieces
59150	Orange	16	180 mm	5.3 mm	Female	30 pieces
59151	Red	18	180 mm	6.0 mm	Female	30 pieces
59152	Green	6	280 mm	2.0 mm	Pediatric	30 pieces
59153	Blue	8	280 mm	2.6 mm	Pediatric	30 pieces
59154	Black	10	280 mm	3.3 mm	Pediatric	30 pieces