




Coated soft tip



USF

For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or for therapeutic purposes.

CHARACTERISTICS

Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining

- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter
- Polished side eyes and tube flexibility additionally reduce friction between the catheter and urethral mucosal lining Soft, flexible tip enables easy catheter introduction in case of enlarged prostate and various urethral masses













SOFT



Catalogue No. Colour Length **Outer diameter** Packaging CH Note 62549 ■ Blue 8 400 mm 2.6 mm 30 pieces Male 62550 ■ Black 30 pieces 10 400 mm 3.3 mm Male 62551 □White 12 400 mm 4.0 mm Male 30 pieces Green 14 4.6 mm 30 pieces 62552 400 mm Male 62553 Orange 16 400 mm 5.3 mm Male 30 pieces 62554 ■ Blue 8 180 mm 2.6 mm Female 30 pieces 62555 3.3 mm ■ Black 10 180 mm Female 30 pieces 62556 □White 12 180 mm 4.0 mm Female 30 pieces 180 mm 4.6 mm 62557 Green 14 Female 30 pieces 62558 Orange 16 180 mm 5.3 mm 30 pieces Female 62560 Blue 8 310 mm 2.6 mm Pediatric 30 pieces 62561 ■ Black 10 310 mm 3.3 mm 30 pieces Pediatric