



Cultura[®] M

A small and handy incubator for
dipslides and contact slides

- Small
- Durable
- Reliable

Cultura[®] M

Cultura M, a small, durable and reliable incubator is suitable for the incubation of Orion Diagnostica's contact slide Hygicult[®] and dip slides Uricult[®], Dentocult[®], and Easicult[®].

The transparent door allows you to check inside the incubator without opening the door. Cultura M is small and light and can be easily transported wherever you need to use it.

Cultura M is equipped with a thermometer and multi-rack shelf. Racks are available also separately.



CE

This CE Mark indicates that this equipment complies with the following directive:

- EMC 89/336/EEC
- DIN EN 61010-4
- 61010-2010

Cultura M Incubator

Temperature range	25...45°C
Temperature stability	± 1°C
Voltage	230 V, 115 V, 12 V
Frequency	50/60 Hz
Protective Class	II
Dimensions (mm)	external W: 310 H: 155 D: 168 internal W: 220 H: 120 D: 150
Case	Flame resistant PC & ABS
Door	PC, transparent
Thermometer	up to 60°C, in special tube
Base cover	ABS
Warranty	2 years from initial installation
Weight	1.1 kg

Ordering information

P/N		Description
05866	129997	Cultura M, 220/230 V
06266	132754	Cultura M, 110/115 V
-	133192	Cultura M, 220/230 V with UK plug
68237	132697	Cultura S, 12 V, with a car-plug
05724	129981	Thermometer, accessory
-	139248	Multirack, accessory