

Gipssäge "HippOscill" 2002 Electronic Power
Plaster Saw "HippOscill" 2002 Electronic Power



| | |
|---|---|
| Art.-No. 1.920.00 | Art.-No. 1.920.01 |
| 230 V~, 50 Hz | 110 V~, 60 Hz |
| Power 500 Watt | Power 500 Watt |
| 6500 – 24000 oscillations, infinitely variable | 6500 – 24000 oscillations, infinitely variable |
| 5 meter cable with EURO-plug | 3 meter cable with USA-plug |
| <p>Plaster saw complete with carton 1 piece sawblade Ø50 mm for normal plaster 1 piece sawblade Ø65 mm for normal plaster 1 pair wrenches for sawblade changing 1 piece instruction manual</p> | |
| <p>Optional (extra charges) : Version in front, however with sawblades for synthetic dressings Convenient, practical, robust and pleasant form case for transport- and / or filing incl. inlay Identification plate service Sawblades in many different forms and coatings for varied employments</p> | |

Gipssäge "HippOscill 2002 Elektronik Power" Plaster Saw "HippOscill 2002 Electronic Power"

- ✓ **500 Watt power**
- ✓ **6500 bis 24000 oscillations, infinitely variable**
- ✓ **Fullautomatical power-compensation above the whole oscillations field**
- ✓ **Slim motor edition with silent motor 52 – 77 dB(A) in relation to the selected oscillations**
- ✓ **Latest electronic system with various protection and security functions**
- ✓ **Special positioning of the motor, protected by plaster dust**
- ✓ **Robust construction for highest loading capacity by continuous operation**
- ✓ **Safety ON / OFF switch**
- ✓ **Weight only 1,25 kg without cable (1,5 kg with cable)**
- ✓ **Vibrations- and maintenance-free**
- ✓ **Regulation from outside of the manual operating area**
- ✓ **5 meter cable with EU-connection, 230 Volt version or 3 meter cable with US-connection, 110 Volt version**
- ✓ **Double protection isolation according to class II, medical use of typ BF**
- ✓ **Radio and TV suppressed**
- ✓ **CE conformity**
- ✓ **EN 60601-1, EN 60601-1-2, IEC 60601-1, IEC 60601-1-2 inspected**

Anton Hipp GmbH

Instruments & Implants

Annastraße 25/1

D-78567 Fridingen

Tel: ++49 7463/99303-0 Fax: ++49 7463/99303-50

E-mail: info@anton-hipp.de www.anton-hipp.de

Distributed by.: